



ASTRUM
INSTRUMENTS



DENTAL
EXTRACTING
FORCEPS



Wazirabad Road Grave Yard Street,
Harrar Near USKT, Slalkot, 51310 Pakistan



+923076616686




Email : astruminstruments@gmail.com
Web : www.astruminstruments.com

EXTRACTING FORCEPS

ENGLISH PATTERN




101
Fig. 1
Upper centrals
& canines




102
Fig. 2
Upper laterals
& bicuspids
& bicuspids




103
Fig. 3
Upper incisors
& canines




104
Fig. 5
Lower incisors
& canines



105
Fig. 7
Upper primolars



106
Fig. 8
Lower primolars



107
Fig. 13
Lower primolars

EXTRACTING FORCEPS

ENGLISH PATTERN



EXTRACTING FORCEPS


ENGLISH PATTERN



115
FIG. 23
RIGHT LOWER MOLARS




116
Fig.24
Left lower molars




117
Fig.29
Upper root & incisors




118
Fig. 29N
Roots of upper incisors



119
Fig.30
Upper roots & incisors




120
Fig.31
Lower roots




121
Fig.32
Lower molars

EXTRACTING FORCEPS


ENGLISH PATTERN




122
Fig.33
Lower roots




123
Fig.33A
Lower roots




124
Fig.33C
Lower roots




125
Fig.33M
Lower roots with long beak




126
Fig.41
Upper roots



127
Fig.44
Upper roots



128
Fig.45
Lower roots




129
Fig.51
Upper roots

EXTRACTING FORCEPS


ENGLISH PATTERN




130
Fig.51A
upper roots narrow beak



131
Fig.51C
Upper roots



132
Fig.52
Upper roots & canines
with wider beak




133
Fig.52A
Upper teeth & resection


115
FIG. 23
RIGHT LOWER MOLARS




134
Fig.55
Separating upper molars roots



135
Fig.56
Separating lower molars roots




136
Fig.67
Upper wisdoms




137
Fig.67A
Upper wisdoms

EXTRACTING FORCEPS


ENGLISH PATTERN




138
Fig. 68
Lower roots




139
Fig. 73
Lower molars (hawk's bill)




140
Fig. 74
Lower incisors & roots



141
Fig. 74M
Lower incisors & roots



142
Fig. 74N
Lower incisors & roots



143
Fig. 75
Lower premolars & canines




144
Fig. 76
Upper roots




145
Fig. 76N
Upper roots
Narrow break

EXTRACTING FORCEPS

ENGLISH PATTERN



146
Fig.79
Lower wisdoms



147
Fig.79A
Lower wisdoms




148
Fig.86A
Lower molars




149
Fig. 86B
Lower molars




150
Fig.86C
Lower molars



151
Fig.87
Cowhorn beak




152
Fig.89
Right upper molars




153
Fig.9
Left upper molars

EXTRACTING FORCEPS


ENGLISH PATTERN



154
Fig.94
Right upper molars




155
Fig.95
Left upper molars




156
Fig.99
Lower molars
(brpken-down)




157
Fig.99 1/4
Lower molars
(brpken-down)




158
Fig.121
Upper wisdoms



159
Fig.136
Upper centrals
& roots




160
Fig.137
Lower centrals
& roots



161
Fig.167
Universal patter
for lower teeth
molars & premolars


EXTRACTING FORCEPS

ENGLISH PATTERN




EXTRACTING FORCEPS


CHILDREN'S PATTERN




170
Fig. 158
Children's upper teeth molars universal




171
Fig. 159
Children's lower molars universal



172
Fig. 160
Children's upper lower universal (hawk's bill)



173
Fig. 161
Children's lower & roots




174
Fig. 162
Children's upper centrals



175
Fig. 163
Children's upper centrals

EXTRACTING FORCEPS

ENGLISH PATTERN



176

Fig.37


Children's upper incisors & canines



177

Fig.39


Children's upper molars & universal



178

Fig.123

Children's lower teeth universal



179

Fig.138

Children's upper universal & roots



180

Fig.157

Children's upper molars universal

EXTRACTING FORCEPS


ROOT FRAGMENT FORCEPS




Wisdom Teeth Extracting Forceps
FOR CHILDREN'S WITH SPRING




185
Fig. 3
Upper molars




186
Fig. 3
Lower incisors




187
Fig. 5
Lower molars




188
Fig. 7
Lower roots



189
Fig.137
Upper incisors
& canines




190
Fig.139
Upper premolars



191
Fig.51S
Upper roots

Tooth Extracting FORCEPS

AMERICAN PATTERN



192
Fig. 1
Upper centrals, laterals & cuspids

193
Fig. 5
Lower third molars

194
Fig. 6
Splitting forceps


195
Fig. 10 H
Upper molars hook-handle

196
Fig. 10 S
Upper molars straight - handle

197
Fig. 15
Lower molar


Tooth Extracting FORCEPS

ENGLISH PATTERN




Tooth Extracting FORCEPS

AMERICAN PATTERN



204
Fig. 32
Upper molars & pre molars
strong beak




205
Fig. 53 L
Upper molars left side




206
Fig. 53 R
Upper molars right side



207
Fig. 62
Universal forceps
for upper teeth & roots



208
Fig. 65
Upper Incisors & root




209
Fig. 69
Fragments & small roots
Upper & lower

Tooth Extracting FORCEPS


AMERICAN PATTERN



210
Fig. 88 R
Right upper molars




211
Fig. 88L
Left upper molars




212
Fig. 99C
Upper Premolars,
Molar & incisors




213
Fig. 101
Universal forceps
for premolars



214
Fig. 150
Upper premolars,
incisors & roots




215
Fig. 150A
Upper premolars incisors & roots




216
Fig. 150S
Upper premolars incisors & roots
for children

Tooth Extracting FORCEPS


AMERICAN PATTERN




217
Fig. 151
Lower premolars, incisors & roots




218
Fig. 151A
Lower premolars, incisors & roots




219
Fig. 151S
Lower premolars, incisors & roots
for children's



220
Fig. 210H
Upper wisdoms with
hook handle



222
Fig. 222
Lower wisdoms



223
Fig. 286
Upper roots, premolars & incisors